

Fluorescein Sodium Ophthalmic Strips 1 mg, 100/Box

User Guidelines

Consumable

Intended use

Sterile single-use strips impregnated with fluorescein sodium for staining the ocular surface and tear film during eye examination under cobalt blue light.

Intended users

Ophthalmologists, optometrists and trained eye-care staff. For use by trained healthcare professionals. Follow your facility's protocols and applicable regulations.

Read this manual completely before first use and keep it with the product. This document provides general operating, care and safety guidance. Where the product is supplied with manufacturer labeling or instructions, those take precedence. For questions, a replacement copy, or a digital copy of this manual, email sales@scimedstore.com.

1. Safety

Safety warnings and cautions

- For use by trained healthcare professionals. Follow your facility's protocols and applicable regulations.
- Single use only. Discard after one application.
- Do not use if the wrapper is opened or damaged or after the expiry date.
- Remove soft contact lenses before use; fluorescein can stain them.
- Avoid touching the cornea with the strip.
- Use is according to clinical protocol and the practitioner's judgment.

2. Package contents

- 100 individually wrapped fluorescein sodium 1 mg strips

Check contents on receipt. Report missing or damaged items to sales@scimedstore.com within 48 hours of delivery.

3. Specifications

Active ingredient	Fluorescein sodium USP
Amount	1 mg per strip
Pack size	100 strips per box
Sterility	Individually wrapped, single use
Storage	Cool, dry place, away from direct sunlight (15-25 °C per listing)
Shelf life	36 months from manufacture (per listing); see box for expiry

4. Operation

Use

- 1 Check the expiry date and open the wrapper at the handle end.
- 2 Moisten the strip tip with a drop of sterile saline.
- 3 Touch the moistened tip gently to the lower conjunctiva or fornix.
- 4 Allow the dye to spread with blinking.
- 5 Examine with cobalt blue light or the slit lamp blue filter.
- 6 Discard the strip.

5. Cleaning and disinfection

- Strips are single use and are not cleaned or reprocessed.

6. Maintenance and storage

- Store in original box in a cool, dry place away from direct sunlight.
- Do not refrigerate or freeze.

- Use strips in order of expiry.

7. Disposal

- Dispose of used strips as clinical waste per local rules.
- Dispose of expired, unopened strips per local pharmaceutical or clinical waste rules.

8. Warranty and support

Products are covered by SciMed's standard warranty (12 months from purchase unless a different period is stated on your invoice). Warranty terms: scimedstore.com/pages/warranty.

Support, spare parts and additional copies of this manual: sales@scimedstore.com. Please include your order number and a photo of the product or issue.

SciMed Scientific Equipment LLC · San Jose, California, USA · scimedstore.com