

# Ocular Prosthetic Eyes, Pack of 10, Brown Shades

## Instructions for Use

Ophthalmic Prosthetic

### Intended use

Stock acrylic ocular prostheses for clinical inventory, used by ocularists and eye care facilities for trial fitting and dispensing to patients.

### Intended users

For use by trained healthcare professionals (ocularists, ophthalmologists and their clinical staff). Follow your facility's protocols and applicable regulations. Patients should handle the prosthesis only as instructed by their eye care professional.

Read this manual completely before first use and keep it with the product. This document provides general operating, care and safety guidance. Where the product is supplied with manufacturer labeling or instructions, those take precedence. For questions, a replacement copy, or a digital copy of this manual, email [sales@scimedstore.com](mailto:sales@scimedstore.com).

## 1. Safety

### Safety warnings and cautions

- For use by trained healthcare professionals. Follow your facility's protocols and applicable regulations.
- Fitting, modification and suitability for each patient must be determined by a qualified ocularist or physician.
- Check the expiry date and case seal before use. Do not use if the case is damaged or opened, or past the expiry date.
- Inspect each prosthesis for cracks, chips, rough edges or scratches before fitting. Do not fit a damaged prosthesis.
- Wash hands or wear gloves before handling. Handle by the edges over a soft surface to prevent drops.
- Advise patients to contact their eye care professional if they have pain, discharge, redness or swelling of the socket.
- Do not use alcohol or solvents on the prosthesis; they can damage acrylic.

## 2. Package contents

- Ocular prostheses x 10 in individual cases

Check contents on receipt. Report missing or damaged items to [sales@scimedstore.com](mailto:sales@scimedstore.com) within 48 hours of delivery.

## 3. Specifications

|            |                                                              |
|------------|--------------------------------------------------------------|
| Quantity   | 10 prostheses                                                |
| Material   | Medical-grade acrylic (PMMA)                                 |
| Sizes      | 20, 22, 24, 26, 28 mm (mixed)                                |
| Shades     | Brown shade family (5 tones)                                 |
| Laterality | Mixed left and right                                         |
| Finish     | Hand-painted iris with veining detail                        |
| Packaging  | Individual sterilized cases; expiry date marked on each unit |

## 4. Operation

### Selection and fitting

- 1 Assess the patient's socket and select the size, side and shade according to clinical judgment.
- 2 Check the expiry date and case seal, then open the case.
- 3 Inspect the prosthesis surface and edges.
- 4 Fit and assess the prosthesis according to your clinical protocol.
- 5 Record the size, side, shade and lot or expiry information in the patient record.

### Patient instructions

- 1 Show the patient how to insert and remove the prosthesis as directed by the fitting clinician.
- 2 Explain how to clean the prosthesis and how often, according to your clinical advice.
- 3 Advise the patient to store it in its case when out of the eye.
- 4 Arrange follow-up and routine checks as recommended by the clinician.

## 5. Cleaning and disinfection

- Clean the prosthesis as instructed by the fitting clinician, usually with clean water and a mild, non-abrasive soap, then rinse well.
- Do not use alcohol, solvents, abrasive cleaners or very hot water.
- Do not autoclave or heat the prosthesis.
- Dry with a soft, lint-free cloth.

## 6. Maintenance and storage

- Store unopened units in their sealed cases at room temperature, away from direct sunlight and heat.
- Use stock in order of expiry date.
- Advise patients to have the prosthesis checked and polished as recommended by their ocularist.
- Replace a prosthesis that is scratched, cracked or no longer fits comfortably.

## 7. Disposal

- Dispose of expired, damaged or used prostheses in accordance with facility procedures and local regulations.

## 8. Warranty and support

Products are covered by SciMed's standard warranty (12 months from purchase unless a different period is stated on your invoice). Warranty terms: [scimedstore.com/pages/warranty](https://scimedstore.com/pages/warranty).

Support, spare parts and additional copies of this manual: [sales@scimedstore.com](mailto:sales@scimedstore.com). Please include your order number and a photo of the product or issue.

SciMed Scientific Equipment LLC · San Jose, California, USA · [scimedstore.com](https://scimedstore.com)