

Stainless Steel Kidney Tray Set, 4 Pieces

Instructions for Use

Surgical Instrument

Intended use

Reusable kidney-shaped stainless steel trays for holding instruments, dressings and supplies and for collecting fluids, emesis or wound drainage in clinical, surgical, dental and veterinary settings.

Intended users

Healthcare, dental and veterinary staff. For use by trained healthcare professionals. Follow your facility's protocols and applicable regulations.

Read this manual completely before first use and keep it with the product. This document provides general operating, care and safety guidance. Where the product is supplied with manufacturer labeling or instructions, those take precedence. For questions, a replacement copy, or a digital copy of this manual, email sales@scimedstore.com.

1. Safety

Safety warnings and cautions

- For use by trained healthcare professionals. Follow your facility's protocols and applicable regulations.
- Treat used trays and their contents as potentially infectious.
- Wear gloves when handling used trays.
- Trays become hot during and after autoclaving; allow to cool before handling.
- Do not use trays that are cracked or have sharp, damaged edges.
- Do not leave sharps loose in trays; use a sharps container.

2. Package contents

- 4 kidney trays: 6 in (150 mm), 8 in (200 mm), 10 in (250 mm), 12 in (300 mm)

Check contents on receipt. Report missing or damaged items to sales@scimedstore.com within 48 hours of delivery.

3. Specifications

Sizes	150 / 200 / 250 / 300 mm (6 / 8 / 10 / 12 in)
Material	Medical-grade stainless steel (SS304 typical per listing)
Construction	One-piece pressed
Cover	None (open tray)
Sterilization	Autoclavable per listing

4. Unpacking and setup

- 1 Clean and sterilize trays before first use.
- 2 Choose the tray size suitable for the task.

5. Operation

Use

- 1 Place the sterile or clean tray on a stable surface or hold it close to the patient as needed.
- 2 Use the curved inner edge against the body for collecting fluids.
- 3 Place instruments or supplies in the tray without overfilling.
- 4 After use, carry the tray level to avoid spills.
- 5 Empty contents into the appropriate waste or disposal system.

6. Cleaning and disinfection

- Rinse off blood and debris promptly after use.
- Wash with a neutral detergent and soft brush; avoid steel wool and chlorine-based cleaners.
- Rinse and dry thoroughly.

- Sterilize per facility protocol; the listing states the trays are autoclavable.

7. Maintenance and storage

- Store clean and dry, stacked by size.
- Inspect for stains, corrosion or damage before use.
- Remove stains with a cleaner suitable for stainless steel.

8. Troubleshooting

Problem	Likely cause	What to do
Water spots or stains	Hard water or residue	Rinse with demineralized water and dry fully
Rust spots	Contact with other corroding items or chlorine	Clean with stainless steel cleaner; avoid chlorine products
Tray hot after sterilizing	Normal	Allow to cool before handling

9. Disposal

- Decontaminate before disposal.
- Recycle as stainless steel scrap at end of life.

10. Warranty and support

Products are covered by SciMed's standard warranty (12 months from purchase unless a different period is stated on your invoice). Warranty terms: scimedstore.com/pages/warranty.

Support, spare parts and additional copies of this manual: sales@scimedstore.com. Please include your order number and a photo of the product or issue.

SciMed Scientific Equipment LLC · San Jose, California, USA · scimedstore.com