

10 mL Universal Pipette Tips, Polypropylene

User Guidelines

Consumable

Intended use

Clear polypropylene pipette tips for use with 10 mL pipettors for transferring and dispensing liquids in laboratories.

Intended users

Laboratory staff and researchers.

Read this manual completely before first use and keep it with the product. This document provides general operating, care and safety guidance. Where the product is supplied with manufacturer labeling or instructions, those take precedence. For questions, a replacement copy, or a digital copy of this manual, email sales@scimedstore.com.

1. Safety

Safety warnings and cautions

- Tips are supplied non-sterile; sterilize before use if your work requires sterility.
- Use single use only to avoid cross-contamination.
- Wear gloves and eye protection when handling hazardous liquids.
- Check fit on your pipettor before use to avoid leaks.

2. Package contents

- 100 tips per bag (no rack)

Check contents on receipt. Report missing or damaged items to sales@scimedstore.com within 48 hours of delivery.

3. Specifications

Volume	10 mL (10,000 µL)
Material	Polypropylene, clear
Sterility	Non-sterile; autoclavable
Pyrogenicity	Non-pyrogenic
Packaging	100 tips per bag
Fit	Universal cone, most 10 mL pipettors

4. Operation

Use

- 1 Load tips into a rack or take one directly from the bag with gloved hands.
- 2 Press the pipettor firmly into the tip to seal.
- 3 Aspirate and dispense following your pipettor's instructions.
- 4 Eject the tip into the appropriate waste container.

5. Cleaning and disinfection

- Tips are single use and are not cleaned.
- Autoclave unused tips before use if sterility is required.

6. Maintenance and storage

- Store in a clean, dry place at room temperature.
- Keep the bag closed to keep dust out.
- Keep away from direct sunlight and heat.

7. Disposal

- Dispose of used tips according to the hazard of the liquid handled.
- Autoclave or decontaminate tips used with biological material before disposal per local rules.

8. Warranty and support

Products are covered by SciMed's standard warranty (12 months from purchase unless a different period is stated on your invoice). Warranty terms: scimedstore.com/pages/warranty.

Support, spare parts and additional copies of this manual: sales@scimedstore.com. Please include your order number and a photo of the product or issue.

SciMed Scientific Equipment LLC · San Jose, California, USA · scimedstore.com